

 erylğit







ERYİĞİT, with its powerful staff and expertise is an export oriented manufacturer and vendor of a wide variety of products ranging from ICU beds to lockers; with an emphasis on health sector with 'Medical Equipment and Furniture' activities.

As company headquarters are located in Ankara in order to ensure effective communications with the clientele.

ERYİĞİT offers top quality, endurance, functionality, fast delivery and professional service as well as flexible prices to fit international competitiveness.

In addition to CE Certification in all product lines, **ERYİĞİT** is also awarded ISO 9001:2000 quality management system certificate and TSE (Turkish Standards Institution) certificates with its concerns for total quality.

As research and development activities gain pace to increase product variety, service quality and customer satisfaction, **ERYİĞİT** with stout determination, takes strong steps towards becoming an internationally acclaimed company for its quality products and services, its technology and professionalism.

4 MOTORS ICU BED

1. Legrest and backrest adjustments are by linear actuator, height, trendelenburg and reverse-trendelenburg adjustments by lifting columns.
2. Electronic and Manual CPR.
3. Fowler, Vascular, Cardiac, Legrest Reverse V positions.
4. Ergonomic hand control unit.
5. Retractable nurse control panel at the end of the bed.
6. Locking positions from the nurse control panel.
7. Sliding backrest.
8. 24 V chargable battery.
9. ABS plastic easy movable head and foot boards.
10. ABS plastic movable, retractable and separate siderails.
11. ABS plastic covered mattress base with 4 sections.
12. ABS plastic covered bottom chassis.
13. Height adjustable IV Pole.
14. Drainage hangers on two sides of the bed.
15. Swivel castors with central, total and directional lock.
16. Four side bumpers.
17. Electrostatic powder paint metal parts.

SPECIFICATIONS

LENGTH	222 CM
WIDTH	100 CM
HEIGHT	45-80 CM
BACKREST ANGLE	75°
FOOTREST ANGLE	30°
TRENDELENBURG	18°
REV - TRENDELENBURG	18°
CASTORS Ø	125-150 MM

± 5 cm tolerance



4 MOTORS ICU BED

1. Legrest, backrest, height, trendelenburg and reverse-trendelenburg adjustments by linear actuator.
2. Electronic and Manual CPR.
3. Fowler, Vascular, Cardiac, Legrest Reverse V positions.
4. Ergonomic hand control unit.
5. Retractable nurse control panel at the end of the bed.
6. Locking positions from the nurse control panel.
7. Sliding backrest.
8. 24 V chargable battery.
9. ABS plastic easy movable head and foot boards.
10. ABS plastic movable, retractable and separate siderails.
11. ABS plastic covered mattress base with 4 sections.
12. ABS plastic covered bottom chassis.
13. Height adjustable IV Pole.
14. Drainage hangers on two sides of the bed.
15. Swivel castors with central, total and directional lock.
16. Four side bumpers.
17. X-Ray conductive headrest.
18. Electrostatic powder paint metal parts.

SPECIFICATIONS

LENGTH	222 CM
WIDTH	100 CM
HEIGHT	45-80 CM
BACKREST ANGLE	75°
FOOTREST ANGLE	30°
TRENDELENBURG	18°
REV - TRENDELENBURG	18°
CASTORS Ø	125-150 MM

± 5 cm tolerance



4 MOTORS ICU BED

1. Legrest, backrest, height, trendelenburg and reverse-trendelenburg adjustments by linear actuator.
2. Electronic and Manual CPR.
3. Fowler, Vascular, Cardiac, Legrest Reverse V positions.
4. Ergonomic hand control unit.
5. Retractable nurse control panel at the end of the bed.
6. Locking positions from the nurse control panel.
7. Sliding backrest.
8. 24 V chargable battery.
9. ABS plastic easy movable head and foot boards.
10. ABS plastic movable, retractable and separate siderails.
11. ABS plastic covered mattress base with 4 sections.
12. ABS plastic covered bottom chassis.
13. Height adjustable IV Pole.
14. Drainage hangers on two sides of the bed.
15. Swivel castors with central, total and directional lock.
16. Four side bumpers.
17. Electrostatic powder paint metal parts.

SPECIFICATIONS

LENGTH	222 CM
WIDTH	100 CM
HEIGHT	45-80 CM
BACKREST ANGLE	75°
FOOTREST ANGLE	30°
TRENDELENBURG	18°
REV - TRENDELENBURG	18°
CASTORS Ø	125 MM

± 5 cm tolerance



SPECIFICATIONS

LENGTH	222 CM
WIDTH	100 CM
HEIGHT	45-80 CM
BACKREST ANGLE	75°
FOOTREST ANGLE	30°
TRENDELENBURG	18°
REV - TRENDELENBURG	18°
CASTORS Ø	125 MM

± 5 cm tolerance

4 MOTORS ICU BED

1. Legrest, backrest, height, trendelenburg and reverse-trendelenburg adjustments by linear actuator.
2. Manual CPR.
3. Fowler, Vascular, Cardiac, Legrest Reverse V positions.
4. Ergonomic hand control unit.
5. Retractable nurse control panel at the end of the bed.
6. Locking positions from the nurse control panel.
7. Sliding backrest.
8. 24 V chargable battery.
9. ABS plastic easy movable head and foot boards.
10. ABS plastic movable, retractable and separate siderails.
11. ABS plastic covered mattress base with 4 sections.
12. ABS plastic covered bottom chassis.
13. Height adjustable IV Pole.
14. Drainage hangers on two sides of the bed.
15. Swivel castors with diagonal lock.
16. Four side bumpers.
17. Electrostatic powder paint metal parts.



4 MOTORS ICU BED

1. Legrest, backrest, height, trendelenburg and reverse-trendelenburg adjustments by linear actuator.
2. Fowler, Vascular, Cardiac, Legrest Reverse V positions.
3. Ergonomic hand control unit.
4. Sliding backrest.
5. ABS plastic easy movable head and foot boards.
6. ABS plastic movable, retractable and separate siderails.
7. ABS plastic covered mattress base with 4 sections.
8. ABS plastic covered bottom chassis.
9. Height adjustable IV Pole.
10. Drainage hangers on two sides of the bed.
11. Swivel castors with diagonal lock.
12. Four side bumpers.
13. Electrostatic powder paint metal parts.

SPECIFICATIONS	
LENGTH	222 CM
WIDTH	100 CM
HEIGHT	45-80 CM
BACKREST ANGLE	75°
FOOTREST ANGLE	30°
TRENDELENBURG	18°
REV - TRENDELENBURG	18°
CASTORS Ø	125 MM

± 5 cm tolerance



SPECIFICATIONS

LENGTH	222 CM
WIDTH	100 CM
HEIGHT	45-80 CM
BACKREST ANGLE	75°
FOOTREST ANGLE	30°
CASTORS Ø	125 MM

± 5 cm tolerance

3 MOTORS HOSPITAL BED

1. Legrest, backrest and height adjustments by linear actuator.
2. Fowler, Vascular and Legrest Reverse V positions.
3. Ergonomic hand control unit.
4. Sliding backrest.
5. ABS plastic easy movable head and foot boards.
6. Movable, retractable and separate steel siderails.
7. ABS plastic covered mattress base with 4 sections.
8. ABS plastic covered bottom chassis.
9. Height adjustable IV Pole.
10. Drainage hangers on two sides of the bed.
11. Swivel castors with diagonal lock.
12. Four side bumpers.
13. Electrostatic powder paint metal parts.



SPECIFICATIONS

LENGTH	222 CM
WIDTH	100 CM
HEIGHT	45-80 CM
BACKREST ANGLE	75°
FOOTREST ANGLE	30°
CASTORS Ø	125 MM

± 5 cm tolerance

3 MOTORS HOSPITAL BED

1. Legrest, backrest and height adjustments by linear actuator.
2. Fowler, Vascular and Legrest Reverse V positions.
3. Ergonomic hand control unit.
4. Sliding backrest.
5. ABS plastic easy movable head and foot boards.
6. ABS plastic movable, retractable and separate siderails.
7. ABS plastic covered mattress base with 4 sections.
8. ABS plastic covered bottom chassis.
9. Height adjustable IV Pole.
10. Drainage hangers on two sides of the bed.
11. Swivel castors with diagonal lock.
12. Four side bumpers.
13. Electrostatic powder paint metal parts.



SPECIFICATIONS

LENGTH	222 CM
WIDTH	100 CM
HEIGHT	50 CM
BACKREST ANGLE	75°
FOOTREST ANGLE	30°
CASTORS Ø	125 MM

± 5 cm tolerance

2 MOTORS HOSPITAL BED

1. Legrest and backrest adjustments by linear actuator.
2. Fowler, Vascular and Legrest Reverse V positions.
3. Ergonomic hand control unit.
4. Sliding backrest.
5. ABS plastic easy movable head and foot boards.
6. ABS plastic movable, retractable and separate siderails.
7. ABS plastic covered mattress base with 4 sections.
8. Height adjustable IV Pole.
9. Drainage hangers on two sides of the bed.
10. Swivel castors with diagonal lock.
11. Four side bumpers.
12. Electrostatic powder paint metal parts.



2 MOTORS HOSPITAL BED

1. Legrest and backrest adjustments by linear actuator.
2. Fowler, Vascular and Legrest Reverse V positions.
3. Ergonomic hand control unit.
4. Sliding backrest.
5. ABS plastic easy movable head and foot boards.
6. ABS plastic movable, retractable and separate siderails.
7. ABS plastic covered mattress base with 4 sections.
8. ABS plastic covered bottom chassis.
9. Height adjustable IV Pole.
10. Drainage hangers on two sides of the bed.
11. Swivel castors with diagonal lock.
12. Four side bumpers.
13. Electrostatic powder paint metal parts.

SPECIFICATIONS

LENGTH	222 CM
WIDTH	100 CM
HEIGHT	50 CM
BACKREST ANGLE	75°
FOOTREST ANGLE	30°
CASTORS Ø	125 MM

± 5 cm tolerance



THREE ADJUSTMENTS MANUAL PATIENT BED

1. Backrest adjustment between the angle of 00-75° and legrest adjustment between the angle of 00-30° by lever.
2. Trendelenburg position with gas spring.
3. Fowler, Vascular, Cardiac, Legrest Reverse V positions.
4. ABS plastic easy movable head and foot boards.
5. ABS plastic movable, retractable and separate siderails. separate siderails.
6. ABS plastic covered mattress base with 4 sections.
7. Height adjustable IV Pole.
8. 4 castors with diagonal brake.
9. Four side bumpers.
10. Drainage hangers on two sides of the bed.
11. Electrostatic powder paint metal parts.
12. Dimensions : L222xW100xH50 (+/- 5cm)



OPERATING ROOM TRANSFER STRETCHER

1. Consist of 2 stretchers and 1 top couch. 2. The couch's backrest can be manually adjustable between the angle of 0°-70°. 3. Stretchers have Trendelenburg position with gas spring. 4. Stretchers can be locked and detached easily. 5. Oxygen tube holder on the stretchers. 6. Tray for materials on the bottom. 7. Four castors with diagonal brake. 8. Four side bumpers. 9. Height adjustable IV Pole. 10. Movable, retractable and lockable CrNi siderails. 11. Artificial leather upholstery covered on 50 mm foam can be detachable for easy cleaning. 12. Pushing handles. 13. Drainage hangers on two sides of the stretchers. 14. Electrostatic powder paint metal parts. 15. Dimensions : L423xW77xH80 (+/- 5cm)



OPERATING ROOM TRANSFER STRETCHER WITH

1. Consist of 2 stretcher and 1 top couch. 2. The couch's backrest can be manually adjustable between the angle of 0°-70°. 3. Stretchers have Trendelenburg position with gas spring. 4. Stretchers can be locked and detached easily. 5. Oxygen tube holder on the stretchers. 6. Tray for materials on the bottom. 7. Swivel castors with central, total and directional lock. 8. Four side bumpers. 9. Height adjustable IV Pole. 10. Movable, retractable and lockable CrNi siderails. 11. Artificial leather upholstery covered on 50 mm foam can be detachable for easy cleaning. 12. Pushing handles. 13. Drainage hangers on two sides of the stretchers. 14. Electrostatic powder paint metal parts. 15. Dimensions : L423xW77xH80 (+/- 5cm)



TKF 900



PATIENT STRETCHER

1. The backrest can be manually adjustable between the angle of 0°-70°.
2. Trendelenburg position with gas spring.
3. Oxygen tube holder on the stretcher.
4. Tray for materials on the bottom.
5. Pushing handles.
6. Movable, retractable and lockable CrNi siderails.
7. Artificial leather upholstery covered on 50 mm foam can be detachable for easy cleaning.
8. Height adjustable IV Pole.
9. Drainage hangers on two sides of the stretcher.
10. Four castors with diagonal brake.
11. Four side bumpers.
12. Electrostatic powder paint metal parts.
13. Dimensions : L214xW71xH65 (+/- 5cm)

TKF 920



EMERGENCY STRETCHER (HYDRAULIC)

1. The backrest can be manually adjustable between the angle of 0°-70°.
2. Height adjustment by hydraulic pump.
3. Foot pedal on both sides.
4. Trendelenburg and Rev-Trendelenburg positions.
5. Four castors with diagonal brake.
6. Pushing handles.
7. Artificial leather upholstery covered on 50 mm foam can be detachable for easy cleaning.
8. Movable, retractable and lockable CrNi siderails.
9. Drainage hangers on two sides of the stretcher.
10. Four side bumpers.
11. Height adjustable IV Pole.
12. Electrostatic powder paint metal parts.
13. Dimensions : L214xW71xH55-80 (+/- 5cm)

TKF 700

EXAMINATION COUCH

1. The backrest can be adjustable between the angle of 0°-70°.
2. Towel holder at the head side.
3. Slots for IV stands both at the head and foot sides.
4. Artificial leather upholstery covered on 50 mm foam.
5. Raising of the head section.
6. With IV holder.
7. Electrostatic powder paint metal parts.
8. Dimensions : L185xW60xH65 (+/- 5cm) (without upholstery)

TKF 710

EXAMINATION COUCH (WITH BEDSTAND)

1. The backrest can be adjustable between the angle of 0°-70°.
2. One bedstand under the couch with castors.
3. Towel holder at the head side.
4. Slots for IV stands both at the head and foot sides.
5. Artificial leather upholstery covered on 50 mm foam.
6. Raising of the head section.
7. With IV holder.
8. Electrostatic powder paint metal parts.
9. Dimensions : L185xW60xH65 (+/- 5cm) (without upholstery)

**TKF 720**

EXAMINATION COUCH (WITH 2 BEDSTANDS)

1. The backrest can be adjustable between the angle of 0°-70°.
2. Two bedstands under the couch with castors.
3. Towel holder at the head side.
4. Slots for IV stands both at the head and foot sides.
5. Artificial leather upholstery covered on 50 mm foam.
6. Raising of the head section.
7. With IV holder.
8. Electrostatic powder paint metal parts.
9. Dimensions : L185xW60xH65 (+/- 5cm) (without upholstery)



TKF ST 100 - TKF SA 100

ST 100 STOOL WITH BACKREST

- 1- Height adjustable
- 2- Plastic base.
- 3- 5 Antistatic swivel castors.
- 4- Artificial leather upholstery.
- 5- Dimensions: L57xW57xH102-115 cm (± 10 cm)



SA 100 IV HANGER

- 1- Height adjustable with two hangers.
- 2- Plastic base.
- 3- 5 Antistatic swivel castors.
- 4- CrNi Steel frame.
- 5- Dimensions: L57xW57xH120-210 cm (± 10 cm)



TKF PA 100 - 110

FOLDING SCREEN (DOUBLE)

- 1- Have two wings.
- 2- Foldable.
- 4- The curtains are detachable for cleaning and they are washable.
- 5- Electrostatic powder coating steel frame.
- 6- Dimensions: W140xH171 cm (± 10 cm)



FOLDING SCREEN (TRIPLE)

- 1- Have three wings.
- 2- Foldable.
- 4- The curtains are detachable for cleaning and they are washable.
- 5- Electrostatic powder coating steel frame.
- 6- Dimensions: W210xH171 cm (± 10 cm)



LOCKERS



TKF-SD 100 *
190x45x42 cm
Single



TKF-SD 110 *
190x45x42 cm
Single & Drawer



TKF-SD 120 *
190x85x42 cm
2 Doors



TKF-SD 130 *
190x85x42 cm
2 Doors & Drawer



TKF-SD 140 *
190x85x42 cm
4 Doors

Bedstands and Dining Tables

TKF-KM 120 *
Dimensions: 75x45x40 cm.



TKF-KM 130 *
Dimensions: 75x45x40 cm.



TKF-KM 140 *
Dimensions: 75x45x40 cm.



TKF-KM 110 *
Dimensions: 75x45x40 cm.



TKF-KM 150 *
Dimensions: 75x45x40 cm.



TKF-KM 100 *
Dimensions: 75x45x40 cm.



TKF-YM 100 *
Top: 40x80 cm.
Height: 83-107 cm



TKF-YM 110 *
ABS Plastic Top : 40x82 cm.
Height: 83-107 cm



TKF-ES 100 *
Dimensions: 45x38x40 cm.

VISCO ELASTIC MATTRESS

TKF HY 150

Dimensions
190x90 cm







QUALITY CERTIFICATES

- ISO 9001:2000
 - CE
 - TSE
- TSE SERVICE PLACE CONFIRMATION
 - QUALITY CERTIFICATE
- AFTER SALES SERVICE CERTIFICATE
 - INDUSTRY REGISTRATION

Our Factory has the right to change
the products and their specifications
without any notice.



www.eryigit.com.tr
foreigntrade@eryigit.com.tr



Ankara - Factory :

İvogsan Özanadolu San. Sit. 1453. Sokak. No.3
TR 06370 Ostim / Yenimahalle / Ankara / Türkiye
P: +90 312 395 57 95 (pbx) - Faks: +90 312 395 57 96



ISO-9001:2000
Quality Management System